

Product datasheet

Human Beta2 Glycoprotein 1 Antibody IgM ELISA Kit

Catalog # EK00832E

Overview

Product name Human Beta2 Glycoprotein 1 Antibody IgM ELISA Kit

Size 96 wells

Species Reactivity Human

Assay Type Sandwich

Sample Type Serum, Plasma, Other human body fluids & secretions

Assay Range $25\mu g/L - 850\mu g/L$

Kit Components Assay Plate, Standard, Standard diluent, Sample diluent, Enzyme

conjugate, Detection reagent A & B, Wash buffer, Stop solution,

Protocol, Quality control certificate

Storage Store at 2-8 °C

Product Note This ELISA Kit uses the double-antigen sandwich method to determine

the level of human anti- β 2 glycoprotein 1 antibody IgM (β 2-GP1 IgM) in the specimen. The ELISA analytical biochemical technique of the kit is

based on β2 GP1 antigen-β2 GP1 IgM antibody interactions

(immunosorbency) and an HRP colorimetric detection system to detect $\beta 2\text{-}GP1$ antibody IgM targets in samples. $\beta 2\text{-}GP1$ IgM are added to the antigen-coated microwells, and then combine with the HRP-labeled antigen to form an antigen-antibody -Enzyme-labeled antigen complex, after thorough washing, add substrate TMB to develop color. TMB is converted into blue under the catalysis of HRP enzyme, and into the final

yellow under the action of acid. The color intensity is positively correlated with $\beta 2$ glycoprotein 1 antibody IgM ($\beta 2$ -GP1 IgM) in the sample. Absorbance (OD value) is measured at a wavelength of 450nm

with a microplate reader, and the concentration of human $\beta 2$ glycoprotein 1 antibody IgM ($\beta 2$ -GP1 IgM) in the sample can be

calculated from the standard curve.







TEL:1-661-524-0262



Typical Data

Results of a typical standard curve are provided with the product for demonstration only and should not be used to obtain test results.

A standard curve must be run for each set of samples assayed.

IMPORTANCE

After receiving the kit, please check the intactness of outer packaging, check whether the reagent bottles are broken, leaking or the amount of liquid is significantly reduced, etc. If the diluent is polluted such as flocculent bacterial colonies; the substrate is not transparent and the microtiter plate is leaking, etc; please contact the distributor. In rare cases, crystals may appear in the substrate. Normally the crystals will dissolve at room temperature, this will not affect the testing results.





TEL:1-661-524-0262